



ALL-PRO 1400 SPRAY

NEUTRAL SURFACE DISINFECTANT

- ALCOHOL & QUAT BASED - READY TO USE

DESCRIPTION

ALL-PRO 1400 SPRAY is a wide spectrum disinfectant and is efficient against: virus, bacteria, yeast, mold and algae. **ALL-PRO 1400 SPRAY** is approved by **Agriculture Canada** as a ready to use wide spectrum disinfectant. **ALL-PRO 1400 SPRAY** belongs to a family of quaternary compounds that are environmentally friendly and solubilised in a solution of isopropyl alcohol.

DILUTION

ALL-PRO 1400 SPRAY is used as is, no dilution in order to control 99.99% of tested pathogenic microorganisms in 10 minutes. However, **ALL-PRO 1400 SPRAY** is efficient in less 15 seconds.

TEMPERATURE

The variation of temperature does not affect the efficiency of **ALL-PRO 1400 SPRAY**.

ADVANTAGES

ALL-PRO 1400 SPRAY is a neutral disinfectant that prevents corrosion and damage to equipments and surfaces. **ALL-PRO 1400 SPRAY** is not irritating to skin.

APPLICATION

Dental clinics: **ALL-PRO 1400 SPRAY** is used to disinfect medical instruments and surfaces. **ALL-PRO 1400** is used undiluted.